


PHOSPHAT CEMENT

normalhärtend / normal setting


MADE IN GERMANY

REF 80022- Farbe)

REF 8002300

REF 80025- Farbe)

Rev. 2024-02

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Hoffmann Dental Manufaktur GmbH
Wangenheimstraße 37/39
D - 14193 Berlin, Germany
☎: +49(0)30/82 00 990
☎: +49(0)30/82 28 145
quality@hoffmann-dental.com
🌐 hoffmann-dental.com

HOFFMANN'S CEMENT
normalhärtend
Zinkphosphat-Cement

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HOFFMANN'S CEMENT
normal setting
Zinc phosphate cement

Instruction for use for the dentist

Composition:

Powder contains: zinc oxide, magnesium oxide; Liquid contains: o-phosphoric acid.

Areas of Application:

Cavity lining as basis for all filling materials (amalgam, composites etc.), cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as metal ceramics (zircon, aluminium and lithium disilicate ceramics) and crowns and bridges on implants, cementation of orthodontic bands, retention pins and screws, cement build up, temporary fillings.

Dispensing Instructions:

Cementation consistency: 1.5g powder : 1.0g liquid
Cavity lining consistency: 2.1g powder : 1.0g liquid
Dispensing and testing climate for testing in accordance with DIN EN ISO 9917: Mix quantities indicated at 23±1°C and a relative humidity of 50±10%, he dispensing instructions also apply when mixing the liquid with **HOFFMANN'S UNIVERSAL CEMENT** .

Mixing:

A clean, dry glass plate at ca. 20°C (68°F) and a non-discolouring spatula should be used for mixing. We recommend Hoffmann's glass mixing slab and Hoffmann's spatula . In order to achieve a homogenous consistency, divide the measured quantity of powder into 4 portions (1/2 1/4 1/8 1/8). Beginning with the smallest quantity, stroke the powder into the measured liquid within 90 seconds. The cavity lining consistency has been reached when on lifting the spatula, the peak drawn out forms a small hook and does not sink back into the paste (Fig. 1). The cementation consistency has been reached when the peak formed on lifting the spatula slowly falls back into the paste (Fig. 2 and 3). Discard remaining powder or impure liquid.

Mixing Time: 90 s

Processing time from beginning of mixing (at 23°C, 72°F):

For cementation: 3:00 min

For cavity linings: 2:30 min

Setting Time from end of mixing:

For cementation: 5:00 – 7:30 min

For cavity linings: 3:30 – 5:30 min

Special Information:

- It is vital that opened pulp is covered. According to the indication we recommend Hoffmann's PULPINE , Hoffmann's PULPINE NE , Hoffmann's PULPINE MINERAL or Hoffmann's ENDO REPAIR as capping material.
- By deep cavity short-term irritation of the pulp could occur, however, this irritation will recede completely after a few minutes. For cavity lining we recommend to use Hoffmann's Copal Varnish .
- To obtain an optimal pH-value, comply with the mixing time of 90 s.
- The cement should be applied to a preferably dry tooth surface.
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.
- Close containers immediately after removing material.

Package Sizes:

powder 100g | liquid 40ml | 30g powder + 15ml liquid

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