HOFFMANN'S

READY2PROTECT Copperioncement

CE 0482

Made in Germany

10 x  je 1g / 1x  5 ml

REF 8007010

Stand der Information / Date of issue
06/2026

Hergestellt von/manufactured by

HOFFMANN DENTAL MANUFAKTUR GMBH

Wangenheimstraße 37/39
14193 Berlin
Germany

Tel: +49(0)30/82 00 99 0
Fax: +49(0)30/82 28 14 5
quality@hoffmann-dental.com
www.hoffmann-dental.com

HOFFMANN'S READY2PROTECT

de

Copperioncement

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en

Copperciment

**Bactericidal protection cement with Copperions for luting and as dentin replacement in single doses for individual consistency control, x-ray visible
In accordance with DIN EN ISO 9917**

Instruction for use for the dentist

Areas of Application:

Minimal invasive treatment of caries with ART Technique; Cavity lining with bactericidal effect on biofilm for all types of teeth (also vital tooth) and as basis for all filling materials (amalgam, composites, glass ionomer, etc.), cementation of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as metal ceramics (zircon-, aluminium and lithium disilicate ceramics) and crowns and bridges on implants, cementation of orthodontic bands, retention pins and screws, cement build up, temporary fillings with bactericidal effect especially suitable for milk teeth. Healing cement to preserve vital teeth, also for deep caries in combination with Hoffmann's Copal Varnish.

Dispensing, Consistency and Times:

To comply READY2PROTECT with individual requirements to the consistency, the user has the ability to vary the amount of drops (see table). It is conform to the requirements of the norm DIN EN ISO 9917.

The setting times for the indicated ratio were determined according to DIN EN ISO 9917: Mixing at 23 ± 1 °C, $50 \pm 10\%$ relative humidity and setting at 37 ± 1 °C, relative humidity at least 90%. In practice, however, there might be slight variations to these times, which are, for example, due to different temperature and humidity conditions.

Mixing:

We recommend Hoffmann's mixing block (Ord-Nr. 82900) and Hoffmann's spatula (Ord-Nr. 82902). Apply the required amount of drops of READY2PROTECT liquid onto the glass slab. To do so hold the tube in an upright position. Build up pressure with index finger and thumb and keep upright up to the last drop (attention: no pumping action). Open READY2PROTECT powder stickpack at the perforation and empty completely onto a clean and dry glass slab. For best possible mixing results, we recommend that the powder be divided into single portions (1/4, 1/4, 1/2). The entire powder is mixed into the liquid with a spatula as follows: 2 x 1/4 in 15 seconds each and 1/2 in 30 seconds. The spatula should be moved in a sweeping fashion – not a stirring movement.

Note:

View an instructional film on mixing at
www.hoffmann-dental.com

Special Information:

- To obtain an optimal pH-value, comply with the mixing time of 60s.
- Store at room temperature
- Do not shake the tube before usage
- It is vital that dentin near the pulp or opened pulp is covered.
- By deep cavity short-term irritation of the pulp could occur, however, this irritation will recede completely after a few minutes. For cavity lining we recommend to use Hoffmann's Copal Varnish (Ord-Nr. 81901).
- The cement should be applied to a preferably dry tooth surface
- The liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.
- Do not use after the expiration date

Application	Consistency	Dispensing: Drops of liquid : 1 Stick		Times		
				Mixing Time	Processing Time	Setting Time
Cementation: (enough material up to 3 crowns)		min	6	60 s	< 3:00 min	< 5:00 – 7:30 min
	creamy	recommended	7		> 3:00 min	5:00 – 7:30 min
	dropwise	max	8		3:00 min	> 5:00 – 7:30 min
Cavity lining	packable	recommended	4	60 s	2:30 min	3:30 – 5:30 min
	pasty	max	5		> 2:30 min	> 3:30 – 5:30 min

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