READY2PROTECT

Copperioncement + Kopallack
Copperioncement + Copal Varnish
Copperioncement + vernis de copal

Made in Germany

REF 8007020

Stand der Information / Date of issue
 06/2026

Hergestellt von/manufactured by
HOFFMANN DENTAL MANUFAKTUR GMBH

Wangenheimstraße 37/39
 D - 14193 Berlin
 Germany
 Tel: +49(0)30/82 00 99 0
 Fax: +49(0)30/82 28 14 5
 quality@hoffmann-dental.com
 www.hoffmann-dental.com

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Flüssigkeit: 1 x Tube à 5 ml

Kopallack: 1 Flasche à 5 ml

Healing cement for preservation of vital teeth also for deep caries in combination with Hoffmann's Copal Varnish, x-ray visible, in accordance with DIN EN ISO 9917

Areas of Application:

Minimal invasive treatment of caries with modified ART Technique with copper sealer; cavity lining with bactericidal effect for all types of teeth (also vital tooth) and as basis for all filling materials (composites, amalgams, etc.); luting of inlays, onlays, crowns and bridges made of precious metal, non-precious metal as well as all-ceramics (zircon, aluminum and lithium disilicate ceramics); permanent cementation of crowns and bridges on implants; cementation of orthodontic bands; retention pins and screws; cement build up; temporary fillings with bactericidal effect; milk teeth fillings.

Dispensing, Consistency and Times:

Consistency control of READY2PROTECT is realized by precise liquid dosing. A whole stick pack with one gram of powder is always used. The amount of drops is varied depending on the required consistency i.e. luting and lining (see below). Requirements of DIN EN ISO 9917 are always met.

Setting times for the indicated ratio were determined according to DIN EN ISO 9917: Mixing at $23 \pm 1^\circ\text{C}$, $50 \pm 10\%$ relative humidity and setting at $37 \pm 1^\circ\text{C}$, relative humidity at least 90%. In practice, however, these vary due to, for example, different temperature and humidity conditions.

Mixing:

We recommend Hoffmann's glass mixing slab (REF 82900) and Hoffmann's spatula (REF 82902). Apply required amount of drops of READY2PROTECT liquid onto the glass slab. To do so hold the tube in an upright position. Buildup pressure with index finger and thumb and keep upright up to the last drop (attention: no pumping). Open READY2PROTECT powder stick at perforation and empty completely onto mixing slab. For best possible mixing results, we recommend to divide powder into 3 single portions (1/4, 1/4, 1/2). The entire powder is mixed into liquid with a spatula as follows: 2 x 1/4 in 15 seconds each and 1/2 in 30 seconds. Spatula should be moved in a sweeping fashion – not a stirring movement.

Note:

View an instructional film on mixing at
www.hoffmann-dental.com

Recommendation for the use on vital teeth, in ART Technique, in pediatric and geriatric dentistry:

It is necessary to protect dentin near the pulp with a copper sealer before cementation.
The application of a copper sealer on a vital tooth stump is also recommended before taking impression.

Copper sealer with copal varnish

0.08 g of copal varnish with 0.1 g of READY2PROTECT Copperioncement powder
Mix together powder and copal varnish by using a single-use brush. Apply a thin layer of the mixture on the cavity walls under the use of the brush and let it dry with compressed air. After complete drying of the first layer the application of a second layer can be repeated in the same way.

Special Information:

- Do not store in fridge.
- Do not shake tube before use.
- Liquid contains phosphoric acid and is corrosive. In the event of accidental eye contact, rinse immediately with large quantities of water.
- Hoffmann´s Copal varnish should not be used in the event of known allergies to its components. Avoid contact with skin and mucous membrane.
- Keep containers tightly close and seal immediately after use.
- Avoid contact of copal varnish bottle with open flame.
- Do not use after expiration date.

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